

EKKO PORTABLE

Speech Therapy Device for Children

🏠 Home Therapy Solution 📱 Mobile App Controlled 📶 Wireless Connectivity

i Introduction

EKKO-portable is an innovative therapeutic device specifically designed for speech therapy of children at home. This medical-grade device is completely controlled through a comprehensive mobile application that assists parents in performing effective therapy sessions. The device has been rigorously tested against safety standards **EN 60601-1: 2006+AMD1: 2013** by 3rd Party CTI (Centre Testing International) labs, Shenzhen, China, and the design is approved by Pakistan Innovation and Testing Centre (PITC) HQ, Islamabad.

⚙️ Technical Specifications

Operating Frequency:	Min: 12Hz±2 Max: 95Hz±2
Operating Voltage:	4.2 V
Battery:	Single cell, 3500 mAH
Average Current:	270mA
Dimensions:	5cm × 20cm × 3.835cm (W × L × H)
Wireless Technology:	Bluetooth Low Energy (BLE)
Weight:	Lightweight & Portable
Certification:	EN 60601- 1: 2006+AMD1: 2013

🔧 Hardware Features

- ✓ **Mobile App Controlled:** Complete device operation via smartphone
- ✓ **Bluetooth Interface:** Seamless BLE connectivity with mobile app
- ✓ **Long Battery Life:** 3500mAH for extended therapy sessions
- ✓ **Compact Design:** Portable and child-friendly form factor
- ✓ **Safety Certified:** Medical device standards compliant
- ✓ **Therapeutic Range:** 12Hz-95Hz frequency spectrum



EKKO Portable Device

Mobile Application Features

 **User Authentication:** Secure login system for patient data protection

 **Device Control:** Complete remote control of EKKO device via smartphone

 **Cloud Data Logging:** Automatic therapy session recording and progress tracking

 **Patient Configuration:** Customizable settings based on specific speech disorders

 **Progress Analytics:** Real-time therapy progress monitoring and reports

 **Session Management:** Automated timing and frequency control during therapy

Usage Guide & Safety Information

▶ Getting Started

Step 1: Device Setup

1. Charge the device fully before first use
2. Ensure device is powered on (check power indicator)
3. Keep device within 10 meters of smartphone

Step 2: App Installation

1. Download EKKO app from Google Play or App Store
2. Create account and complete authentication
3. Enable Bluetooth on your smartphone
4. Pair device via Bluetooth in app settings

Step 3: Patient Configuration

1. Input child's age and speech therapy requirements
2. Configure frequency settings as per therapist recommendations
3. Set therapy duration and session schedule
4. Complete initial device calibration through app

⚙ Operating Instructions

Mobile App Operation:

- All device functions controlled via smartphone app
- Real-time frequency adjustments during sessions
- Automated therapy programs based on child's needs
- Cloud synchronization for progress tracking
- Remote monitoring and control capabilities

App Interface Features:

- Current frequency settings display
- Battery level monitoring
- Bluetooth connection status indicator
- Active therapy session timer
- Progress tracking dashboard
- Therapy history and analytics

Session Guidelines:

- Ensure stable Bluetooth connection before starting
- Follow app-guided therapy protocols
- Begin with recommended frequency settings
- Keep initial sessions short (5-10 minutes)
- Monitor child's comfort through app feedback

Important Safety Information:

- ! Always supervise children during therapy sessions
- ! Do not exceed recommended frequency limits (12-95Hz)
- ! Discontinue use immediately if child shows discomfort
- ! Keep device away from water and moisture
- ! Charge only with provided charger

- ! Clean device with soft, dry cloth only

Medical Considerations:

- 🔗 Consult speech therapist before first use
- 🔗 Follow prescribed therapy guidelines strictly
- 🔗 Monitor child's progress regularly via app
- 🔗 Report any adverse reactions immediately
- 🔗 Maintain regular follow-ups with healthcare provider

🔧 Maintenance

Regular Maintenance:

- Clean device weekly with soft, dry cloth
- Check battery level regularly through app
- Update mobile app when notifications appear
- Store in dry, cool environment when not in use
- Avoid exposing device to extreme temperatures

Battery Care:

- Charge when battery level shows below 20% (Red LED ON)
- Avoid complete battery discharge
- Use original charger and cable only
- Expected battery life: 2-3 years with normal use

Common Issues & Solutions:

Device won't power on:

- Check battery charge level
- Try hard reset (contact support for instructions)
- Ensure proper charger connection

Bluetooth connection fails:

- Restart both device and smartphone
- Clear Bluetooth cache in phone settings
- Re-pair device through EKKO app
- Ensure devices are within 10-meter range

Mobile app issues:

- Update app to latest version
- Restart smartphone application
- Check internet connection for cloud features
- Clear app cache if problems persist

For Technical Support & Inquiries

✉ contact@ekko.pk 📞 +92 343 9309980 🌐 <https://ekko.pk/ekko-portable/>

Certified by Pakistan Innovation and Testing Centre (PITC) |
EN 60601-1: 2006+AMD1: 2013 Compliant